CERAFIX HV

● NOTICE D'INSTRUCTIONS

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English

Associated instructions:

Reference	Name
646	High viscosity cement without antibiotic

Before using Cerafix HV, the surgeon should familiarise him/herself with the following points:

1. RECOMMENDATIONS

- To ensure safe and efficient use of CERAFIX® bone cements, the surgeon must receive specific training in their use and be familiar with their handling, usage limitations and application technique.
- As each component is pre-dosed, the whole of the solid component should be mixed with the whole of the liquid component.
- No product should be added or placed in contact during the preparation and application of this cement, to avoid any polymerisation problems or any alteration of its physical and mechanical characteristics.
- The use of iodine derivatives or rifamycin or rifampicin-class antibiotics may slow or even inhibit the polymerisation reaction.
- As the liquid monomer is volatile and flammable, the operating theatre must be correctly ventilated and all precautions related to the use of flammable products must be taken.
- Any contact between the monomer and the skin, mucosal surfaces or latex objects (surgical gloves) must be avoided.
- Prolonged inhalation of monomer vapours may have a soporific effect. Excessive exposure to concentrated vapours may also irritate the respiratory tract and eyes.
- Concentrated monomer vapours may have adverse effects on soft contact lenses.


2. DESCRIPTION

The cement is obtained by mixing two pre-dosed sterile single-use components, one a liquid and the other a solid.


Liquid component

With a volume of 21 ml contained in an amber glass ampoule, it has the following composition:

- | | |
|-------------------------------|---------|
| • Methyl methacrylate | 17.7 ml |
| • N-butyl methacrylate | 2.9 ml |
| • N, N-diméthyl paratoluidine | 0.4 ml |
| • Hydroquinone | 45 ppm |

methylmethacrylate


n-butyl methacrylate

Solid component

Weighing 50g and contained in a polyethylene bag, it possesses the following composition:

- | | |
|--------------------------------|--------|
| • Methyl polymethacrylate | 45.0 g |
| • Zirconium dioxide | 4.6 g |
| • Activatable benzoyl peroxide | 0.4 g |

When the solid and liquid components are mixed together the N, N-dimethyl paratoluidine breaks down the benzoyl peroxide to form free radicals that enable monomer polymerisation. This leads to the binding of the polymer beads. Co-polymerisation of methyl methacrylate with a small amount of n-butyl methacrylate improves the cement's impact resistance. The zirconium dioxide renders the cement radio-opaque.

3. INDICATIONS

CERAFIX HV is a high-viscosity, radio-opaque, resin-based acrylic cement specifically designed for fixing endoprostheses to bones during total or partial arthroplasty of the hip, knee, or any other joint.

4. EXPECTED PERFORMANCE

CERAFIX® HV bone cements are of the self-hardening acrylic type. They are used for placing and fixing metal or polymer prostheses in living bone, acting by mechanically bonding the bone surface to the prosthesis. CERAFIX® HV cements provide stable and durable anchoring of osteoarticular prostheses through an even distribution of stress transmitted to the cement towards the bone when the patient moves. Moreover, the anchoring can be clearly identified on X-rays thanks to the radio-opaque properties of the zirconium dioxide present in the cement.

5. PREPARATION - INSTRUCTIONS FOR USE

- Preparation is performed in a sterile area.
- The cement components must be allowed to reach the ambient temperature of the operating theatre.
- Prepare the following materials and keep them at operating theatre temperature: a dry and sterile cup (polyethylene, ceramic or glass), several sterile spatulas, a cement syringe, if required to inject the cement, and a timer.

MIXING PHASE

- Pour all of the solid component into the cup and mix thoroughly with a spatula.
- Pour all of the liquid component onto the powder. Start the timer.
- Blend the liquid into the solid mass using a spatula.
- Mix carefully until a homogeneous texture is obtained, avoiding air-bubble formation.

WAITING PHASE

- Allow to stand away from all heat sources

CEMENT-USE PHASE (MANUAL USE)


- Take the paste without kneading it as soon as it no longer sticks to the gloves. It should be matt in colour.
- Apply the cement.
- Position the implant
- Remove any excess cement

HARDENING PHASE

- Hold the implant perfectly still.

6. CHRONOLOGY AND CLINICAL USE

Bone cements are heat-sensitive. Any rise and/or fall in temperature (ambient temperature and/or temperature of the cement components and instruments used for mixing it) will affect the cement's characteristics, handling and setting time. Humidity variations may also alter its handling characteristics and setting time.

7. CONTRAINDICATIONS

- If an infection is present, the use of cement with no added antibiotic must be carefully considered and should generally be avoided.
- Allergies to product components.
- Infectious arthritis and inactive infection of the joints or joints to replace, or a history of such infections.
- The use of CERAFIX HV cement is not recommended for vertebroplasty or kyphoplasty procedures.

8. ADVERSE EFFECTS

The following complications have been observed during the use of methyl-methacrylate-based surgical cements.

- pulmonary embolism,
- greater trochanteric bursitis,
- vein thrombosis,
- heterotopic osteogenesis,
- haemorrhage and haematoma,
- heart rhythm disorders and hypotension,
- strokes,
- reduced or poor range of joint movement,
- fracture.

9. WARNINGS: USE DURING PREGNANCY AND BREASTFEEDING

The innocuousness of the bone cement during pregnancy or breastfeeding has not been established. As a precautionary measure, bone cement should not be used during the first three months of pregnancy and should be avoided during the rest of pregnancy, except in the case of life-threatening diseases.

10. PRECAUTIONS FOR USE

- During use, the cement components, along with all instruments with which they may enter into contact, must be kept at the operating theatre temperature.
- The following parameters may affect the handling and setting-time characteristics of the cement: specific temperature of the components and instruments used to prepare the cement, temperature and hygrometry of the operating theatre, preparation conditions, storage time.
- When using CERAFIX® HV for endoprosthetic anchoring, all osteotomised surfaces that will be in direct contact with the cement must be washed and suction-cleaned before application. After polymerisation, the joint edges must be carefully cleaned to avoid soiling the joint surfaces with cement particles.
- During and after cement application, the patient's blood pressure, pulse and breathing must be carefully monitored. A drop in blood pressure, accompanied by reduced cardiac activity and an increased breathing rate, may occur during cement application.
- As the cement polymerisation reaction is exothermic, the product reaches temperatures higher than physiological temperature during the hardening phase.

11. PACKAGING AND STORAGE

One box contains:

- one sterile ampoule (21ml) of sterile liquid component,
- one sterile bag containing 50.0g of sterile solid component.

The cement should be stored in a dry and dark place at a temperature below 25°C.

12. STERILISATION

- The liquid component is filter-sterilised. The ampoule, in single packaging, is ethylene-oxide-sterilised.
- The solid component is packaged in a polyethylene bag, which is placed inside a strippable external bag. The whole package is sterilised by gamma radiation (minimum dose 25 kGy).
- Sterility is guaranteed only if the various component envelopes are intact.
- Sterility is guaranteed for 3 years after sterilisation date, solely provided that packaging remains intact.

13. RE-STERILISATION AND REUSE

The re-sterilisation and/or reuse of the components of a unit of CERAFIX® HV cement is strictly prohibited, since this could expose the patient concerned to contamination that might cause infection. Reusing the components that make up CERAVÉR® bone cement would reduce the performance of the device concerned.

14. DISPOSAL AFTER USE

In the event of product return, disposal is managed by CERAVÉR. If products are not returned, disposal is the responsibility of the distributor or health establishment, pursuant to local regulations.

15. STABILITY PERIOD

Do not exceed the use-by date indicated on the packaging.

16. NOTE

CERAFIX® HV high-viscosity orthopaedic cements received their first CE marking in 2002.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

 **CERAVÉR - Les Laboratoires Osteal Médical**
69 rue de la Belle Etoile
CS 54263 - 95957 Roissy CDG Cedex - France
Tel.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99
WWW.CERAVÉR.COM ; ADV@CERAVÉR.COM

Date of last revision of this document: March 2018

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Catalogue reference
	Expiration date		Do not re-use
	Batch code		NB: check precautions for use regarding important safety-related information, such as warnings and precautions that must be taken
	Sterilised by irradiation		Date of manufacture
	Sterilised by exposure to ethylene oxide		Sterilised using aseptic techniques
	CE marking followed by identification of Notified Body (0459, LNE/G-MED)		Inflammable
	Do not store at temperatures above 25°C / 77°F		Keep away from sunlight
	Do not use if package is damaged		

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