HiFit PS

STERILE R

CE 0459

- **PROTHÈSE DE GENOU TRICOMPARTIMENTAIRE CIMENTÉE**
Notice d'instructions
- **CEMENTED TRICOMPARTMENTAL KNEE PROSTHESIS**
Instructions for use
- **TRIKOMPARTIMENTELLE ZEMENTIERTE KNIEPROTHESE**
Anleitung
- **PRÓTESIS TRICOMPARTIMENTAL DE RODILLA CEMENTADA**
Instrucciones de utilización
- **PROTESI DI GINOCCHIO TRICOMPARTIMENTALE CEMENTATA**
Istruzioni d'uso
- **ТРЕХКОМПОНЕНТНЫЙ ПРОТЕЗ КОЛЕННОГО СУСТАВА
С ЦЕМЕНТНОЙ ФИКСАЦИЕЙ**
Инструкция по эксплуатации

STERILE R

0459

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• **CEMENTED TRICOMPARTMENTAL KNEE PROSTHESIS**
Instructions for use

Instructions for use for the following implants:

Cemented HiFit PS femoral component	15501	15502	15503	15504	15505	15506	15507	15508	15509
	15601	15602	15603	15604	15605	15606	15607	15608	15609
HiFit PS tibial plate	15907	15909	15911	15913	15915	15917	15920	15927	15929
	15931	15933	15935	15937	15940	15947	15949	15951	15953
	15955	15957	15960	15967	15969	15971	15973	15975	15977
	15980	15987	15989	15991	15993	15995	15997	16000	16007
	16009	16011	16013	16015	16017	16020	16027	16029	16031
	16033	16035	16037	16040	16047	16049	16051	16053	16055
	16057	16060							
Cemented HiFit PS tibial base	15701	15703	15705	15707	15709	15711			

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not re-use
	Expiration date		Note: the precautions for use must be consulted for all important safety information, such as warnings and precautions
	Batch code		
	Sterilised by radiation		CE mark followed by identification of notified body (0459, LNE/G-MED)
	Catalogue reference		

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• CEMENTED TRICOMPARTMENTAL KNEE PROSTHESIS - Instructions for use

Before using a tricompartmental knee prosthesis from the cemented HiFit PS range, surgeons must familiarise themselves with the following recommendations:

1- GENERAL POINTS

1.1 The manufacturer and supplier of CERAVÉR® total knee prostheses guarantee manufacturing and materials quality. A prosthesis that has already been used must not be implanted again. Instructions for assembling the various components must be respected, as specified in the commercial documentation.

Using CERAVÉR® prosthetic knee components with components of other origins is strictly prohibited.

1.2 Several criteria are key to the success of a procedure involving the use of a cemented HiFit PS tricompartmental total knee prosthesis:

- Appropriate patient selection for the procedure,
- Observance of the operative technique and correct use of the surgical instruments specially designed by CERAVÉR®.
- Application of strict aseptic conditions.

1.3 The surgeon is responsible for any complications that may arise from erroneous indication and/or use of the instruments associated with the operative technique and/or of poor asepsis; neither the manufacturer nor the supplier of the CERAVÉR® products can be held responsible for these complications.

1.4 The surgeon must be aware of all possible problems involved in implanting a knee prosthesis. The surgical procedure for implanting a CERAVÉR® total knee prosthesis can best be learned at a clinic with surgeons who have acquired expertise in the technique.

1.5 The surgeon must inform patients receiving prosthetic knee joints of the fact that their weight and activity levels may affect the durability and performance of the prosthesis.

2- DESCRIPTION OF THE PROSTHESIS

The cemented HiFit PS (posterior stabilised) tricompartmental total knee prosthesis is a sliding total knee prosthesis with moderate constraints comprised of three prosthetic components: the femoral component, the tibial base and the tibial plate. The femoral implant and the tibial implant are separate from each other; they roll and slide on each other by means of a polyethylene insert. This device is indicated solely for cemented use. Furthermore, the prosthesis may be combined with two additional components: the patellar button and the tibial extender. Whether these components are used or not is decided by the surgeon, based on the characteristics of the patient in question.

CERAVÉR® knee prostheses are manufactured using the latest technologies. Their mechanical and functional characteristics can only be guaranteed if surgical technique and the present instructions are strictly observed.

Components

The cemented HiFit PS femoral component is manufactured from cast cobalt-chromium-molybdenum alloy (CoCrMo - ISO Standard 5832-4).

The cemented HiFit PS tibial component is made up of two elements assembled during surgery:

- The tibial base is manufactured from forged, cast Ti-Al6-V4 alloy (ISO Standard 5832-3),
- The HiFit PS tibial plate is manufactured from ultra-high molecular weight polyethylene (UHMWPE) (ISO Standard 5834-2).

3- INDICATIONS

The components of the cemented HiFit PS tricompartmental knee prosthesis are indicated solely for cemented prosthesis implant among patients with bone maturity presenting one of the following arthritic conditions: arthritis, post-trauma arthritis, rheumatoid arthritis and necrosis, causing daily handicap that was not sufficiently improved by medical treatment or surgery, after an observation period ranging from a few weeks to a few months.

In general, the patient must suffer from severe joint involvement that is unlikely to respond to conservative treatment or another surgical strategy.

The aim of total knee arthroplasty is to restore all or some of a patient's mobility and stability, while minimising or even eradicating pain entirely.

4- USE

The following factors may be extremely important to the success of the procedure:

4.1 Selecting the implant

Careful selection of the right implant size is necessary for correct joint function. The degree of success of a total knee prosthesis implant depends on choosing an implant of the right size, shape and model. This total knee prosthesis requires careful placement and should be reserved for a limited level of functional constraint.

4.2 Choosing candidates for a total knee prosthesis

Concerning the selection of total knee prosthesis candidates, the following factors may be extremely important to procedure success:

- **Patient's weight:** in obese patients, excessive weight may increase the load applied to the prosthesis and can compromise the integrity of the cement layer and/or of the medical device.
- **Patient's activity or occupation:** if the patient's activity or profession requires him/her to walk, run, lift heavy loads and/or make muscular efforts, the resulting forces can lead to a failure of the cement, biological anchoring and/or of the medical device.
- **Senility, mental disorders or substance abuse (alcoholism for example):** these disorders may, in particular, cause the patient to ignore certain precautions and necessary limitations related to the use of the medical device and can potentially compromise the implant's integrity or lead to other complications.
- **Certain degenerative diseases:** in certain cases, the evolution of degenerative diseases may have progressed to such an extent at the time of implant placement that they may significantly reduce the medical device's service life. In these cases, the placement of a total knee prosthesis should only be considered a temporary remedy for pain, or as an intermediate procedure.
- **Sensitivity to foreign bodies:** if sensitivity to the material used is suspected, conduct all appropriate tests prior to selecting the prosthesis material and performing implant placement.
- **Infection:** a local infection, whether recent or chronic, may constitute a contraindication to the implantation of a total knee prosthesis. Take all the necessary precautions when selecting patients presenting with a recent or chronic infection.

4.3 Handling the medical device

- The exposed implant joint surfaces must not be placed in contact with any metallic objects or other hard objects that could potentially damage the surface.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorised by the manufacturer must not be implanted.

• CEMENTED TRICOMPARTMENTAL KNEE PROSTHESIS - Instructions for use

4.4 Component assembly and compatibility

4.4.1 General points

- The cemented HiFit PS femoral component can be combined with patellar buttons (Ref. 5090, 5091, 5092, 5093).
- Do not exchange CERAVÉR® knee components with those of other manufacturers.
- Do not exchange components from the HiFit PS range with components from the CERAGYR® range.
- Do not exchange CERAVÉR® knee components with those of other manufacturers, with the exception of a different brand of patellar component which, in revision cases, may be retained provided that its shape is similar to our implant, i.e. a portion of a sphere. Indeed, in some cases, it is preferable to keep the patella rather than attempting to remove it, as this risks weakening the patella and reducing the effectiveness of the extensor apparatus. Only the surgeon can make this decision as he/she is the only person able to analyse the associated risks and benefits.
- The tibial component is made up of two parts that must be engaged by using a rear to front movement during surgery.
 - HiFit PS tibial plates are not compatible with HERMES PS tibial bases and vice versa.
 - The locking mechanism can only be activated once. Do not attempt to reinsert a previously used tibial insert.
- To check compatibility between different components, refer to the table below:

		TIBIAL COMPONENT SIZES (PLATE + BASE)					
COMPATIBILITY		HERMES PS 0/1/2/3 (small stud)	HERMES PS 2/3/4/5 (large stud)	HERMES PS Revision 1/2	HERMES PS Revision 3/4	HiFit PS A/B/C/D (small stud)	HiFit PS C/D/E/F (large stud)
FEMORAL COMPONENT SIZES	HERMES PS 1/1.5/2/2.5	YES	NO	YES	NO	YES	NO
	HERMES PS 3/3.5/4/5	NO	YES	NO	YES	NO	YES
	HERMES PS Revision 1/2	YES	NO	YES	NO	YES	NO
	HERMES PS Revision 3/4	NO	YES	NO	YES	NO	YES
	HiFit PS 1/2/3/4/5	YES	NO	YES	NO	YES	NO
	HiFit PS 6/7/8/9	NO	YES	NO	YES	NO	YES

4.4.2 Adding an extender

Tibial extenders are optional components for HiFit PS prostheses.

The cemented HiFit PS tibial base may be combined with the following extenders and screws:

- Tibial extender (5080, 5081, 5082, 5083, 5084, 5085, 5086, 5087, 5088, 5089);
- Extender screw (3051).

HiFit PS tibial base shafts have a screw cap (Ref. 6297) that must be removed before connecting an extender to the tibial base.

* CERAGYR: mobile-bearing knee prosthesis not approved in Canada

For HiFit PS tibial bases, tibial extenders are assembled using a tapered socket connected to an assembly screw (Ref. 3051):

- Mount the extender on the base without seating it;
- Engage the assembly screw in the base until it is screwed into the extender;
- Seat the extender using a 500 g hammer, tapping it at least 3 times while supported on a surface that is as rigid as possible, avoiding metal-on-metal contact by inserting a dressing folded at least 3 times;
- Finish screwing in.

5- CONTRAINDICATIONS

General contraindications

The placement of a HiFit PS tricompartmental prosthesis is contraindicated in the following cases:

- Patient's physical condition likely to compromise the anchoring or correct positioning of implants (pronounced osteoporosis, systemic and metabolic disorders causing progressive skeletal deterioration, malignant tumours, widespread bone cysts, other joint impairment, severe deformity). Conditions potentially causing excessive joint stress or overload (neurogenic arthritic conditions, major muscle impairment, multiple joint disabilities, fixed deformity of the ipsilateral hip in a prejudicial position);
- Local sepsis, infection, persistent skin lesion;
- Poor skin coverage over the knee joint;
- Known allergy to any of the materials used;
- Repeated falls, major neurological problems, a mental or neuromuscular disorder that could result in unacceptable risk for the patient and lead to postoperative complications.

Warnings (factors that may shorten the implant's service life)

- A history of osteotomy or malunion;
- High Body Mass Index (>30 kg/m²);
- Stress due to intense physical activity (sport, occupation);
- Major functional deficiency of the collateral ligaments.

6- WARNINGS

6.1 Preoperative warnings

- Surgical technique and preoperative planning for the implantation of a tricompartmental knee prosthesis from the HiFit PS range are the basic principles of good surgical practices for the fitting of a total knee prosthesis. In-depth knowledge of the surgical technique is required. Specific surgical tools must be used and the surgeon must be familiar with them. Check the alignment and cutting guide instruments before surgery. Damaged or bent instruments may lead to incorrect implant positioning, thus compromising the success of the procedure. The surgical technique brochure associated with the type of prosthesis is available from CERAVÉR®.
- Use X-ray templates to estimate the size and positioning of the implant, and its alignment with the HKA axis. A sufficient stock of implants of various sizes must be available at the time of surgery, including implants that are smaller and larger than planned. Carefully check the integrity of all packaging and implants before surgery.
- Do not, under any circumstances, reuse a surgical implant. After placement, if an implant is removed, it must be discarded. Though it may appear intact, it may present small defects and signs of internal stress that could compromise its integrity. Only use new implants and do not alter their configuration before use.

6.2 Perioperative warnings

- Choosing the right implant is crucial. Selection must be made in accordance with the appropriate size for each patient and must take relevant anatomical and biomechanical factors into account. These factors include, in particular: the patient's age, level of physical activity, weight, muscle conditions and bone quality.
- If the implant used is not the appropriate size, if the component is not correctly matched to the bone sections or if the stability of the medical device cannot be guaranteed, loosening, dislocation, regression or fracture of the components may occur. It is therefore particularly necessary to ensure that there is a tight bone/cement/prosthesis interface.
- Do not use CERAVÉR® components, instruments and test prostheses with those of another manufacturer. As dimensional compatibility cannot be guaranteed, using components by other manufacturers could result in a poor outcome. With the exception of the patellar component, do not exchange components from the HERMES or HiFit ranges with components from the CERAGYR range (to check HERMES/HiFit matches, see above table).
- Implant the tibial component in such a way as to guarantee suitable bone cortex coverage in all directions (medial-lateral and anterior-posterior).
- It is important to ensure adequate tibial surface preparation and cleaning to achieve correct prosthesis anchoring. Keep bone ablation to the strict minimum to allow implant placement. Excessive bone ablation can lead to mechanical disorders, bone resorption and, consequently, failure of the procedure as a result of a loose or deformed implant. When preparing and placing the tibial components, ensure correct tibial alignment.
- The use of polymethyl methacrylate (PMMA) cement may help fix, support and stabilise certain bone implants, but it cannot either replace the support function of healthy bone tissue or preclude the need for additional support during the consolidation period. When cement is used to secure the implant, ensure that the cement is supported on all parts of the component included in the cement in order to avoid any possible concentration of stresses that could compromise the success of the procedure.
- Before closing the surgical site, clean carefully to remove all bone fragments, ectopic bone tissue, cement, etc. The presence of foreign particles in the metal/plastic joint interface may lead to excessive wear and/or friction. Ectopic bone tissue and/or bone spurs may cause dislocation or painful and restricted mobilisation. Check joint mobility carefully to detect any early contact or instability.

7- PRECAUTIONS**7.1 Preoperative precautions**

- When a total knee prosthesis implant is being considered, particularly in a young and/or active patient, all aspects of the surgery and the implant should be discussed with the patient before the procedure. This consultation should cover the limitations of joint reconstruction, any patient-specific limitations, the possible consequences of these limitations and thus the absolute requirement to follow the doctor's preoperative instructions.
- During the preoperative period, and though rare, it is necessary to consider and eliminate any risks of allergy and other reactions to the implant materials.
- Assessment of the safety and efficacy of using HiFit PS tricompartmental prostheses in bilateral applications is the remit of the surgeon.

7.2 Perioperative precautions

- Correct implant handling is extremely important. Do not allow the device to touch another metallic or hard surface prior to implant placement. Do not use this device if it displays scratches or other damage that could cause defects and stress potentially compromising the success of the implant.
- Do not remove then reinsert previously used polyethylene inserts.

7.3 Postoperative precautions

Postoperative care is important.

- The patient must be informed of the limitations imposed by a HiFit PS tricompartmental knee prosthesis, warned about weight-bearing on the operated limb and given information about permitted levels of activity. Monitor initial support attempts closely.
- Carefully structure initial postoperative care in such a manner as to maintain joint mobility and prevent any dislocation or thromboembolism.
- Structure postoperative treatment, patient care (dressing changes, bed pan use, etc.) and physical activities in a way that will prevent overloading the operated knee. The patient's age and/or bone quality may require the limited load period to be extended.
- All patients who have a CERAVÉR® tricompartmental knee prosthesis fitted should see an orthopaedic surgeon for postoperative follow-up. This examination is important to detect any signs of prosthetic wear or loosening before the appearance of functional manifestations.
- Perform regular X-ray examinations for comparison to X-rays taken after the procedure in order to detect any medium- and long-term progressive position changes or loosening of the implant(s), or any signs of bone retraction, component/cement fracture and/or component disassembly.
- When annual examination is impossible, the surgeon should receive, at least once a year, a control X-ray plate for evaluation and to allow early detection of any lesions.
- The implants are compatible with MRI examination. They may nonetheless create artefacts that should be taken into account for a perfect interpretation of the examination. However, a period of 48 hours is recommended between implant placement and MRI examination to avoid any movement of the device.
- Encourage the patient to notify his/her doctor immediately of any unusual changes to the operated joint.
- HiFit PS implant safety, compatibility, heating and migration have not been assessed by MRI.

8- ADVERSE EFFECTS

- Complications of anaesthesia;
- Cement-induced shock;
- Vascular, muscular or ligament lesions or regressive or definitive nerve problems;
- Perioperative and postoperative fracture;
- Haemorrhage, heterologous or autologous transfusion;
- Misalignment of prosthetic components;
- Other complications associated with general surgical procedures, (digestive and cardiovascular problems), with drugs or ancillary surgical equipment used, with blood, etc.
- Haematoma;
- Thromboembolic complications: phlebitis and pulmonary embolism;
- Slow wound healing;
- Early and late-onset infections;
- Inflammation;
- Algodystrophy;
- Ectopic ossification;
- Persistent knee pain, reduced range of joint movement;
- Instability of knee joint, femoropatellar or femorotibial dislocation or sub-dislocation;
- Allergy to any of the implant components;
- Osteolysis-inducing reaction to products generated by corrosion or wear;
- Mobilisation or detachment of prosthetic components;
- Revision surgery.

9- PACKAGING AND STORAGE

- The prostheses are stored in their original packaging.
- The various components are single-use, they are unit packaged in double bags and provided sterile.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorised by the manufacturer must not be implanted.
- Do not use the device after the expiration date. Contact your CERAVÉR® representative to return the product.

10- STERILISATION

- All CERAVÉR® knee implants are sterilised using gamma radiation at a minimum dose of 25 kGy; UHMWPE components are argon-conditioned to avoid polymer oxidation during sterilisation.
- Before use, it is important to check the integrity of the sterility protection.

11- RESTERILISATION AND REUSE

The resterilisation and/or reuse of any of the components of a CERAVÉR® knee prosthesis is strictly forbidden because of the risk of contamination that could lead to infection. The resterilisation and/or reuse of any of the components of a CERAVÉR® knee prosthesis may compromise the device's biological and mechanical performance.

12- DISPOSAL AFTER USE

If a device is removed for follow-up replacement, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

13- NOTE

The HiFit PS range obtained CE marking in 2012. Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

CERAVÉR® - Les Laboratoires Osteal Médical
 69 rue de la Belle Etoile
 CS 54263 - 95957 Roissy CDG Cedex - France
 Tel.: +33 (0)1 48 63 88 63 - Fax: +33 (0)1 48 63 88 99
WWW.CERAVÉR.COM

Date of last revision of this document: April 2015

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