LES LABORATOIRES OSTEAL MEDICAL

HERMES PS REVISION

STERILE R

CE 0459

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STERILE R

• CEMENTED HERMES PS REVISION TRICOMPARTMENTAL KNEE PROSTHESIS

Instructions for use

Instructions for use for the following components:

Left femoral component		5701	5702	5703	5704			
Right femoral component		5801	5802	5803	5804			
Tibial base plate	Monoblock augment 10mm - Total	5241	5242	5243	5244			
	Monoblock augment 10 mm Left medial / Right lateral	5341	5342	5343	5344			
	Monoblock augment 10 mm Right medial / Left lateral	5441	5442	5443	5444			
Tibial insert	Size 1	1196	1197	1198	1199	1200	4921	1201
	Size 2	1202	1203	1204	1205	1206	4922	1207
	Size 3	1208	1209	1210	1211	1212	4923	1339
	Size 4	1214	1215	1216	1217	1218	4924	1213
	Size 2 Peg size 3	1344	1341	1342	1343	1338	4926	1385
	Size 3 Peg size 2	1345	1346	1347	1348	1384	4927	/
Posterior femoral augment	Thickness 5 mm / Right lateral - Left medial	5541	5542	5543	5544			
	Thickness 10 mm / Right lateral - Left medial	5641	5642	5643	5644			
	Thickness 5 mm / Right medial - Left lateral	5551	5552	5553	5554			
	Thickness 10 mm / Right medial - Left lateral	5651	5652	5653	5654			

• **HERMES PS REVISION** - Instructions for use

Distal femoral Augment	Thickness 10 mm / left medial	5561	5562	5563	5564		
	Thickness 10 mm / left medial	5661	5662	5663	5664		
	Thickness 5 mm / Left lateral	5571	5572	5573	5574		
	Thickness 10 mm / Left lateral	5671	5672	5673	5674		
	Thickness 5 mm / Right medial	5581	5582	5583	5584		
	Thickness 10 mm / Right medial	5681	5682	5683	5684		
	Thickness 5 mm / Right lateral	5591	5592	5593	5594		
	Thickness 10 mm / Right lateral	5691	5692	5693	5694		
Screw for the Femoral Augments		5033					
Femoral extension	Angle 5°	5034	5035	5036	5037	5038	5039
	Angle 7°	5044	5045	5046	5047	5048	5049
Tibial stem offset		5806	5807	5808	5809		

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken
	Batch code		
	Sterilized by irradiation		CE marking followed by the code for the notified body (0459, LNE/G-MED)
	Catalogue reference		

Before using a cemented HERMES PS REVISION tricompartmental knee prosthesis, the attention of the surgeon is drawn to the following recommendations:

1 - GENERAL POINTS

- A. The manufacturer and supplier of a CERAVER total knee prosthesis guarantee the quality of the manufacturing process and materials used. A prosthesis which has already been used must not be implanted again. Instructions for the assembly of the different components have to be respected, as specified in these instructions and the commercial documentation.
The use of CERAVER prosthetic knee components with components of other origins is strictly prohibited.
- B. Several criteria are key to the success of a procedure to implant a HERMES PS REVISION tricompartmental total knee prosthesis:
 - Appropriate patient selection for the procedure;
 - Correct implant selection;
 - Observance of the operative technique and correct use of the surgical instruments specially designed by CERAVER;
 - Application of strict aseptic technique.
- C. The surgeon is responsible for any complications that may arise from erroneous indication and/or use of the instruments associated with the operative technique and/or of poor aseptic practice; neither the manufacturer nor the supplier of the CERAVER products can be blamed for these complications;
- D. The surgeon must be aware of all problems posed by implanting a knee prosthesis. The surgical procedure for implantation of a CERAVER total knee prosthesis can be best learned at a clinic or training centre with surgeons who have acquired expertise in the technique;
- E. The surgeon must inform patients receiving prosthetic knee joints of the fact that their weight and activity levels may influence the durability and performance of the prosthesis.

Important note: All users of the HERMES PS REVISION system acknowledge that they have read and accepted the conditions set out in this document, which are to be taken as contractually binding.

2 - DESCRIPTION OF THE PROSTHESIS

The cemented HERMES PS REVISION tricompartmental total knee prosthesis is a total knee sliding prosthesis with moderate stresses comprising three prosthetic components: the femoral component, the tibial baseplate and tibial insert.

The femoral implant and tibial implant are independent from each other; they roll and slide over each other by means of a polyethylene insert.

This device is exclusively indicated for a cemented use.

In the event of a bone defect, the HERMES PS REVISION range offers:

- a removable solution on the femur. The femoral component can be combined with distal and/or posterior filler augments of 5 and 10 mm thickness, fixed by means of a screw;
- a single block solution on the tibia. The tibial baseplate has a monoblock augment of 10 mm thickness, internal, external or total.

The tibial and femoral extensions provide better bone anchoring of femoral and tibial components, often necessary in revision procedures.

Furthermore, the prosthesis can be associated to an additional component: the patellar component. The use or non-use of this component is left up to the assessment of the surgeon according to the characteristics of the patient.

• HERMES PS REVISION - Instructions for use

Components

The cemented HERMES PS REVISION femoral component comprises at least the first two elements assembled during the procedure:

- A HERMES PS REVISION femoral component, made of cobalt chromium-molybdenum alloy (CoCrMo - ISO 5832-4 standard);
- A femoral extension, made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard);
- One or two distal augment(s) (optional), made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard);
- One or two posterior augment(s) (optional), made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard);
- One or more screws for the femoral augments (optional), made of cobalt chromium-molybdenum alloy (CoCrMo - ISO 5832-12 standard).

The cemented HERMES PS REVISION tibial component comprises at least the first three elements listed below, assembled during the procedure:

- A HERMES PS REVISION tibial extension, made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard);
- The HERMES PS REVISION tibial insert, made from ultra high molecular weight polyethylene (UHMWPE - ISO 5834-2 standard);
- A tibial extension, made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard);
- A tibial stem offset adapter (optional), made of titanium alloy (Ti-Al6-V4 - ISO 5832-3 standard).

3 - INDICATIONS

The components of the cemented HERMES PS REVISION tricompartmental knee prosthesis are only indicated for implanting with cement of prostheses for bone defects limited to the epiphysis (below 20 mm at the femur and 15 mm at the tibia) and/or with no collateral ligament failure or slight to moderate failure:

- in the context of replacing an implant, unicompartmental prosthesis (femorotibial or femoropatellar) or sliding tricompartmental prosthesis;
- as a first-line treatment, in specific clinical situations, in patients with mature bones with one of the following arthropathies: arthrosis, post-traumatic arthrosis, rheumatoid arthritis and necrosis causing daily handicap insufficiently improved by medical or surgical treatment, after an observation period of a few weeks to a few months requiring the use of:
 - on the femoral side: an extension and/or augments;
 - on the tibial side: an extension with or without tibial stem offset and/or augments.

4 - ANTICIPATED PERFORMANCE

- A. The aim of partial arthroplasty of the knee is to allow patients to regain stability and ensure the operated knee has a good level of mobility while minimizing and even eliminating pain;
- B. A minimum IKS score of 70 is anticipated;
- C. On the basis of test according to standard ISO 21536, the HERMES PS REVISION prosthesis can have an amplitude of joint movement of 136° before implantation.

5 - CONTRAINDICATIONS

HERMES PS REVISION tricompartmental prosthesis implantation is contraindicated in the following cases:

General contraindications

- A. A physical condition that might compromise adequate implant positioning and fixation (significant osteoporosis, systemic or metabolic problems leading to progressive bone breakdown, malignant tumors, extensive bone cysts or any other joint weakness or severe malformation). Conditions potentially imposing excessive joint stress or overload (e.g.: neurogenic arthropathy, major muscle or ligament insufficiency, multiple joint disabilities, fixed deformity of the overlying hip in an abnormal posture);
- B. History of falls, major neurological disorders, mental or neuromuscular disorders that may constitute an unacceptable risk for the patient and cause postoperative complications;
- C. Known allergy to any of the materials used;
- D. Septic environment, infection and/or persistent skin lesion;
- E. Poor skin coverage over the knee joint.

Specific contraindications:

- F. Metaphyseal-diaphyseal bone defect;
- G. Significant functional defect of collateral ligaments that may reduce the stability of the prosthesis.
- H. Replacement of a hinged prosthesis.

6 - WARNINGS

6.1 Preoperative warnings

6.1.1 General points

- A. A knee prosthesis should be considered only after all other surgical methods of treatment have been carefully weighed up against each other and none have been judged to be more appropriate.
- B. Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good joint mobility and weight-bearing capacity of the operated limb.
- C. Every artificial joint is always subject to wear which cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may render a re-operation necessary.
- D. Infection of a prosthesis will usually necessitate its removal and any such infection may leave debilitating sequelae.

6.1.2 Warnings

Factors that may shorten the lifetime of the implant:

- A. High Body Mass Index (>30 kg/m²): in obese patients, excessive weight may increase the load applied to the prosthesis and compromise the integrity of the layer of cement and/or the medical device.
- B. High-impact physical activity (sport, work activity, etc.): if patients' occupation or activities require them to walk, run, lift weights and/or put strain on their muscles, the resulting forces may cause a failure of the cement, the biological fixation and/or the medical device.
- C. Major misalignment in the frontal and sagittal plane: due to changes in the bone morphology at the ends and/or the presence of osteosynthesis material.

• **HERMES PS REVISION** - Instructions for use

- D. Vascular insufficiency in the affected limb requiring an opinion from a vascular disease specialist.
- E. History of infection of the affected joint and/or history of systemic infection that may affect the prosthesis.
- F. General drug therapy that may lead to a gradual deterioration of the prosthesis' bone support.

6.2 Perioperative warnings

6.2.1 General points

- A. Careful selection of the right implant size is necessary for correct joint function. The degree of success of a total knee prosthesis implantation depends on choosing an implant of the right size and model. The patient's age, level of physical activity, weight, muscle condition and bone quality must also be considered as well as the condition of the collateral ligaments.
- B. If the implant used is not of the optimal size, if the component is not correctly matched to the bone cuts, or if the stability of the medical device cannot be guaranteed, loosening, dislocation, regression or fracture of the components may occur. It is therefore necessary to ensure, in particular, that there is a tight bone/prosthesis or bone/cement/prosthesis interface.
- C. The exposed implant joint surfaces must not be placed in contact with any metallic objects or other hard objects that could potentially damage the surface.
- D. Joint prostheses must not be subjected to treatment or modification.
- E. Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorized by the manufacturer must not be implanted.
- F. Do not remove and reinsert the polyethylene inserts once they are in place.**

6.2.2 Assembly and compatibility of components

- A. The cemented HERMES PS REVISION femoral component can be combined with cemented patellar components, ref.: 5090, 5091, 5092, 5093.
- B. HERMES PS REVISION tibial extensions are not compatible with the tibial inserts of the HiFit PS range.**
- C. HERMES PS REVISION implants are compatibles with HERMES PS first-line treatment prostheses so as to limit repeat operative work when changing the only part of the prosthesis that requires it (femur or tibia) and thus preserve the patient's bone stock.
- D. Do not use CERAVAR devices with those of other manufacturers, apart from, as part of a revision, the patellar component (and provided that it is spherical in shape). The conservation of the patellar component is left up to the surgeon's assessment.
- E. Refer to the following table for the compatibility between the different components of the HERMES PS and HiFit PS range:

COMPATIBILITY		TIBIAL COMPONENT SIZE (INSERT + BASEPLATE)							
		HERMES PS REVISION Size 1 and 2	HERMES PS REVISION Size 3 peg 2 intermediate	HERMES PS REVISION Size 3 and 4	HERMES PS REVISION Size 2 peg 3 intermediate	HIFIT PS A/B (small peg)	HIFIT PS C/D (small peg)	HIFIT PS C/D (large peg)	HIFIT PS E/F (large peg)
FEMORAL COMPONENT SIZE /	HERMES PS REVISION 1	YES	YES	NO	NO	YES	YES	NO	NO
	HERMES PS REVISION 2	YES	YES	NO	YES	YES	YES	NO	NO
	HERMES PS REVISION 3	NO	NO	YES	YES	NO	NO	YES	YES
	HERMES PS REVISION 4	NO	NO	YES	YES	NO	NO	YES	YES
	HIFIT PS 1/2/3/4/5	YES	YES	NO	YES	YES	YES	NO	NO
	HIFIT PS 6/7/8/9	NO	NO	YES	YES	NO	NO	YES	YES

F. The tibial component comprises two parts that must be assembled by clicking them together from the back to the front during the surgical operation.

- **HERMES PS tibial inserts are not compatible with the tibial base plates of the HIFIT PS ranges and vice versa;**
- **The snap mechanism can only be used once. Do not attempt to reinsert a tibial insert that has already been used.**

G. Combination with a femoral augment

One or more posterior and/or distal augments can be associated with the HERMES PS REVISION femoral component:

- Position the augment on the femoral component;
- Screw in the femoral augment screw (ref. 5033).

H. Association of a femoral extension

A femoral extension must be associated with the HERMES PS REVISION femoral component:

- Assemble the extension on the femoral component;
- Strike the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;

I. Association of a tibial extension

The cemented HERMES PS REVISION tibial baseplate must be associated with the extensions (5080, 5081, 5082, 5083, 5084, 5085, 5086, 5087, 5088, 5089) and the screw for tibial extension.

- Assemble the extension in the baseplate without knocking it;
- Insert the assembly screw into the baseplate until it is screwed into the extension;

• HERMES PS REVISION - Instructions for use

- Strike the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;
- Finish screwing.

An offset adapter (ref. 5806, 5807, 5808, 5809) can be inserted between the extension and the tibial baseplate. The tibial extension is then assembled with the offset adapter by conical sleeve:

- Once the offset adapter is positioned in relation to the tibial baseplate, insert the assembly screw;
- Position the offset adapter, as defined during the operation, in the baseplate without knocking it;
- Insert the assembly screw into the baseplate until it is screwed into the offset adapter;
- Place the extension in the offset adapter;
- Strike the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;
- Finish screwing.

7 - PRECAUTIONS

7.1 Preoperative precautions

- When implanting a total knee prosthesis is being considered, particularly in a young and/or active person, all aspects of the surgery and the implant should be discussed with the patient before the procedure. This consultation should cover the limitations of joint reconstruction, any patient-specific limitations, the possible consequences of these limitations and thus the absolute requirement to follow the doctor's preoperative instructions.
- During the preoperative period, it is necessary to consider and eliminate any risks of allergy and other reactions to the implant materials, however rare these risks may be.
- Surgical technique and preoperative planning for the implantation of a HERMES PS REVISION range tricompartmental knee prosthesis correspond to the basic principles of good surgical practices for fitting a total knee prosthesis. An in-depth knowledge of the surgical technique is required.
- Specific surgical tools must be used and the surgeon must be familiar with them. Cutting alignment and guiding instruments should be checked before surgery. Damaged or bent instruments may lead to incorrect implant positioning, thus compromising the success of the procedure. The technical surgical brochure associated with the type of prosthesis is supplied with each set of instruments and is available from CERAVÉR.
- Radiographic templates should be used to estimate the size and positioning of implants as well as its alignment to the HKA axis. A sufficient stock of implants of various sizes must be available at the time of surgery, including sizes smaller and larger than planned. The integrity of all packaging and implants must be carefully checked before surgery.
- It is up to the surgeon to gauge the safety and efficacy of bilateral HERMES PS REVISION range tricompartmental prosthesis implantation.

7.2 Perioperative precautions

- A. If a contraindication comes to light during the operation, it is up to the surgeon to decide whether or not to continue with implantation or implant another type of prosthesis.
- B. The tibial component should be implanted so as to ensure adequate coverage of the bone's cortical substance in all directions (medial-lateral and anterior-posterior).
- C. It is important to ensure adequate tibial surface preparation and cleaning to achieve correct prosthesis anchoring. Keep bone ablation to the strict minimum to allow implant placement. Excessive bone ablation can lead to mechanical disorders, bone resorption and, consequently, failure of the procedure as a result of loosening or deformation of the implant. When preparing and placing the tibial components, correct tibial alignment should be checked.
- D. The use of polymethyl methacrylate (PMMA) cement can help fix, support and stabilize certain bone implants but it cannot either replace the support function of healthy bone tissue or preclude the need for additional support during the consolidation period. When cement is used, complete support of the cement on all parts of the component included in it should be checked in order to prevent the concentration of stresses that could compromise the success of the procedure.
- E. Before closing the surgical site, the knee should be cleaned carefully to remove all bone fragments, ectopic bone tissue, cement, etc. The presence of foreign particles in the metal/plastic joint interface may lead to excessive wear and/or friction. Ectopic bone tissue and/or bone spurs may cause dislocation or painful and restricted mobilization. Check joint mobility carefully to detect any early contact or instability.

7.3 Postoperative precautions

- A. The patient must be told about the restrictions associated with a HERMES PS REVISION range tricompartmental knee prosthesis and warned about using the operated limb for support and informed of permitted activity levels. Monitor initial support attempts closely.
- B. Carefully structure initial postoperative care in such a manner as to maintain joint mobility and to prevent any dislocation or thromboembolism.
- C. Structure postoperative treatment, patient care (dressing changes, bed pan use, etc.) and physical activities such as to prevent any overload of the operated knee. The patient's age and/or bone quality may require the limited load period to be extended.
- D. All patients implanted with a CERAVAR tricompartmental knee prosthesis should see an orthopaedic surgeon for postoperative follow-up. This examination is important to detect any signs of abnormal prosthetic wear or loosening before the appearance of functional manifestations.
- E. Encourage the patient to notify his/her doctor immediately of any unusual changes to the operated joint.
- F. When annual examination is impossible, the surgeon should receive, at least once a year, a control radiograph for evaluation.
- G. The implants are compatible with MRI examination. They may nonetheless create artefacts that should be taken into account for a perfect interpretation of the examination. However, a period of 48 hours is recommended between implant placement and MRI examination to avoid any movement of the device. Safety, compatibility, heating and migration of HERMES PS REVISION implants have not been tested in an MRI. An assessment has been performed in comparison to similar devices available on the market.

8 - ADVERSE EFFECTS

The most commonly reported potential complications are:

- A. Anesthesia complications;
- B. Bone cement implantation syndrome;
- C. Regressive or permanent vascular, musculo-ligamentous or nerve lesions;
- D. Per and postoperative fractures;
- E. Haemorrhage, heterologous/autologous transfusion;
- F. Misalignment of prosthetic components;
- G. Other complications associated with a general surgical procedure (gastrointestinal, cardiovascular, pulmonary, urological problems, etc.), with the medicines or ancillary surgical equipment used, with blood, etc.
- H. Haematoma, sanguineous or non-sanguineous intra-articular effusions;
- I. Thromboembolic complications: deep vein thrombosis and pulmonary embolism;
- J. Wound healing disorders;
- K. Early and late-onset infections;
- L. Inflammation;
- M. Algodystrophy;
- N. Ectopic ossification;
- O. Persistent knee pain;
- P. Reduced or low amplitude of movement;
- Q. Knee joint instability, femoropatellar or femorotibial dislocation or subluxation;
- R. Allergy to implant components;
- S. Reaction to products of corrosion, free particles of cement or wear that may cause osteolysis;
- T. Movement, loosening or breakage of prosthesis components;
- U. New operation.

9 - PACKAGING AND STORAGE

- A. The prostheses are stored in their original packaging, intact and not opened. No special storage conditions
- B. The various components are single-use, they are unit packaged in double bags and provided sterile. The packaging ensures that the product is sealed and not contaminated, given normal handling and transport conditions.
- C. The integrity of the packaging must be checked before use. Do not use if the packaging is damaged Do not use the device after the expiry date. Contact your CERAVER representative to return the product.

10 - STERILIZATION

- A. All CERAVER knee implants are sterilized using gamma radiation at a minimum dose of 25 kGy; UHMWPE components are argon-conditioned to avoid polymer oxidation during sterilization.
- B. Before use, it is important to check the integrity of the sterility protection.

11 - RESTERILIZATION AND REUSE

- A. Resterilization and/or reuse of any of the components of a CERAVÉR knee replacement is strictly forbidden because of the risk of contamination that could lead to infection.
- B. Resterilization and/or reuse of any of the components of a CERAVÉR knee replacement may compromise the device's mechanical performance and biological safety. Never reuse surgical implants. If an implant is removed after it has been implanted, it must be disposed of. Although it may appear intact, it may exhibit minor faults and signs of internal stresses that may undermine its integrity. Only new implants should be used and their configuration must not be modified before use.

12 - DISPOSAL AFTER USE

If a device is removed for revision surgery, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

13 - COMPLAINTS

Any healthcare professional (for example a surgeon using the products) who has a complaint or who has noted a defect relating to the quality, identification, reliability, safety, effectiveness and/or performance of the HERMES PS REVISION system must inform CERAVÉR, or the authorised distributor. In the event of a serious incident, or risk of a serious incident, with the potential to cause or have caused death or severe harm to the health of a patient or a user, CERAVÉR or the authorised distributor must be informed immediately by telephone, fax or mail. For all complaints, we ask that the name(s), reference(s), and lot number(s) of the implicated components be supplied, as well as the name and address of the complainant, the nature of the complaint with as much detail as possible and an indication of the expected response.

14 - NOTE

The HERMES PS REVISION range was granted CE Marking in 1997.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

CERAVÉR - Les Laboratoires Osteal Médical
69 rue de la Belle Etoile
95957 Roissy CDG Cedex - France
Tel.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99
www.ceraver.com ; adv@ceraver.com

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