# HERMES UNI

STERILE R

CE 0459

- **PROTHÈSE DE GENOU UNICOMPARTIMENTAIRE FÉMORO-TIBIALE CIMENTÉE HERMES UNI**  
Notice d'instructions
- **CEMENTED UNICOMPARTMENTAL TIBIOFEMORAL KNEE PROsthESIS HERMES UNI**  
Instructions For Use
- **ZEMENTIERTE, UNIKOMPARTIMENTELLE FEMOROTIBIAL-KNIEPROTHESE HERMES UNI**  
Anleitung
- **PRÓTESIS DE RODILLA UNICOMPARTIMENTAL FEMOROTIBIAL CEMENTADA HERMES UNI**  
Manual de instrucciones
- **PROTESI DI GINOCCHIO MONOCOMPARTIMENTALE FEMORO-TIBIALE CEMENTATA HERMES UNI**  
Istruzioni per l'uso
- **ОДНОМЫШЦЕЛКОВЫЙ ФЕМОРОТИБИАЛЬНЫЙ ПРОТЕЗ КОЛЕННОГО СУСТАВА ЦЕМЕНТНОЙ ФИКСАЦИИ HERMES UNI**  
Инструкция по эксплуатации

STERILE R

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STERILE R

## • CEMENTED UNICOMPARTMENTAL TIBIOFEMORAL KNEE PROSTHESIS HERMES UNI

### Instructions For Use

Instructions for use for the following components:

<b>HERMES UNI cemented femoral component</b>	1228	1229	1230	1231	1302						
<b>HERMES UNI cemented tibial baseplate</b>	1267	1268	1269	1270	1271	1232	1233	1234	1235	1236	
<b>HERMES UNI tibial insert</b>	1237	1238	1239	1243	1244	1245	1249	1250	1251	1255	
	1256	1257	1261	1262	1263	1272	1273	1274	1278	1279	
	1280	1284	1285	1286	1290	1291	1292	1296	1297	1298	

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken
	Batch code		
	Sterilized by irradiation		CE marking followed by the code for the notified body (0459, LNE/G-MED)
	Catalogue reference		

• **HERMES UNI** - Instructions for use

**Before using a HERMES UNI unicompartmental tibiofemoral knee prosthesis, the attention of the surgeon is drawn to the following recommendations:**

**1 - GENERAL POINTS**

- A. The manufacturer and supplier of a HERMES UNI unicompartmental tibiofemoral knee prosthesis from CERAVÉR guarantee the quality of the manufacturing process and materials used. A prosthesis which has already been used must not be implanted again. Instructions for the assembly of the different components have to be respected, as specified in these instructions and the commercial documentation. The use of CERAVÉR prosthetic knee components with components of other origins is strictly prohibited.
- B. Several criteria are key to the success of a procedure to implant a HERMES UNI unicompartmental tibiofemoral prosthesis:
  - Appropriate patient selection for the procedure;
  - Correct implant selection;
  - Observance of the operative technique and correct use of the surgical instruments specially designed by CERAVÉR;
  - Application of strict aseptic technique.
- C. The surgeon is responsible for any complications that may arise from erroneous indication and/or use of the instruments associated with the operative technique and/or of poor asepsis; neither the manufacturer nor the supplier of the CERAVÉR products can be blamed for these complications.
- D. The surgeon must be aware of all problems posed by implanting a knee prosthesis. The surgical procedure for implantation of a HERMES UNI unicompartmental tibiofemoral knee prosthesis can be best learned at a clinic or training center with surgeons who have acquired expertise in the technique.
- E. **The surgeon must inform patients receiving prosthetic knee joints of the fact that their weight and activity levels may influence the durability and performance of the prosthesis.**

**2 - DESCRIPTION OF THE PROSTHESIS**

The HERMES UNI unicompartmental tibiofemoral knee prosthesis includes three prosthetic components: the femoral component, the tibial insert and the tibial baseplate.

CERAVÉR knee prostheses are manufactured using the latest technologies. The mechanical and functional properties can only be guaranteed if the surgical technique and present instructions are strictly observed.

***Components***

- **The femoral component** of the HERMES UNI unicompartmental tibiofemoral knee prosthesis is made of cast cobalt chromium-molybdenum alloy (CoCrMo - NF ISO 5832-4 standard).
- **The tibial component** of the HERMES UNI unicompartmental tibiofemoral knee prosthesis consists of two components that are assembled during the implantation procedure: a metal baseplate and a fixed polyethylene tibialinsert:
  - The tibial baseplate is made of titanium alloy (Ti6Al4V - NF EN ISO 5832-3 standard);
  - The tibial insert is made of ultra high molecular weight polyethylene (UHMWPE) (NF ISO 5834-2 standard).

### 3 - INDICATIONS

The above-mentioned components of the HERMES UNI unicompartmental tibiofemoral knee prosthesis are designed exclusively for implantation with cement in patients with mature bones with one of the following arthropathies located in one of the femorotibial compartments: primary arthrosis, necrosis and post-traumatic arthrosis causing everyday disability that has not improved sufficiently after an observation period of several weeks to several months.

In general, patients must be suffering from severe joint involvement which is unlikely to respond to conservative treatment or another surgical strategy.

### 4 - ANTICIPATED PERFORMANCE

- A. The aim of partial arthroplasty of the knee is to allow patients to regain stability and ensure the operated knee has a good level of mobility while minimizing and even eliminating pain;
- B. A minimum IKS score of 70 is anticipated;
- C. Since the cruciate ligaments are preserved and the prosthesis has an unconstrained design, the joint should be at least as mobile as it was pre-operation.

### 5 - CONTRAINDICATIONS

HERMES UNI unicompartmental prosthesis implantation is contraindicated in the following cases:

- A. A physical condition that might compromise adequate implant positioning and fixation (significant osteoporosis, systemic or metabolic problems leading to progressive bone breakdown, malignant tumors, extensive bone cysts or any other joint weakness or severe malformation). Conditions potentially imposing excessive joint stress or overload (e.g.: neurogenic arthropathy, major muscle or ligament insufficiency, multiple joint disabilities, fixed deformity of the overlying hip in an abnormal posture);
- B. History of falls, major neurological disorders, mental or neuromuscular disorders that may constitute an unacceptable risk for the patient and cause postoperative complications;
- C. Known allergy to any of the materials used;
- D. Active infection, sepsis and/or persistent skin lesion;
- E. Poor skin coverage over the knee joint.

#### Specific HERMES UNI contraindications:

- F. Inflammatory rheumatism;
- G. Impaired function of the cruciate and/or collateral ligaments;
- H. Unreducible flexion contracture or major recurvatum;
- I. Severe arthritis of the contralateral femorotibial compartment and/or femoropatellar compartment;
- J. Major bone misalignment in the frontal plane.

## **6 - WARNINGS**

### **6.1 Preoperative warnings**

#### **6.1.1 General points**

- A. A knee prosthesis should be considered only after all other surgical methods of treatment have been carefully weighed up against each other and none have been judged to be more appropriate.
- B. The HERMES UNI unicompartmental femorotibial prosthesis must be implanted by orthopedic surgeons who regularly perform knee replacement surgery.
- C. Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity of the operated limb.
- D. Every artificial joint is always subject to wear which cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may render a re-operation necessary.
- E. Infection of a prosthesis will usually necessitate its removal and any such infection may leave debilitating sequelae.

#### **6.1.2 Warnings**

Factors that may shorten the implant's lifetime:

- A. High Body Mass Index (>30 kg/m<sup>2</sup>): in obese patients, excessive weight may increase the load applied to the prosthesis and compromise the integrity of the layer of cement and/or the medical device.
- B. High-impact physical activity (sport, work activity, etc.): if patients' occupation or activities require them to walk, run, lift weights and/or put strain on their muscles, the resulting forces may cause a failure of the cement, the biological fixation and/or the medical device.
- C. History of osteotomy or malunion: due to changes in the bone morphology at the ends and/or the presence of osteosynthesis material.
- D. Vascular insufficiency in the affected limb requiring an opinion from a vascular disease specialist.
- E. History of infection of the affected joint and/or history of systemic infection that may affect the prosthesis.
- F. General drug therapy that may lead to a gradual deterioration of the prosthesis' bone support.
- G. Chondrocalcinosis.
- H. Replacing one unicompartmental prosthesis with another is not recommended.

## 6.2 Preoperative warnings

### 6.2.1 General points

- A. Careful selection of the right implant size is necessary for correct joint function. The degree of success of a unicompartamental knee prosthesis implantation depends on choosing an implant of the right size and model. The patient's age, level of physical activity, weight, muscle condition and bone quality must also be considered as well as the condition of the cruciate and collateral ligaments.
- B. If the implant used is not of the optimal size, if the component is not correctly matched to the bone cuts, or if the stability of the medical device cannot be guaranteed, loosening, dislocation, regression or fracture of the components may occur. It is therefore necessary to ensure, in particular, that there is a tight bone/prosthesis or bone/cement/prosthesis interface.
- C. The exposed implant joint surfaces must not be placed in contact with any metallic objects or other hard objects that could potentially damage the surface.
- D. Joint prostheses must not be subjected to treatment or modification.
- E. Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorized by the manufacturer must not be implanted.

### 6.2.2 Assembly and compatibility of components

- A. The cemented HERMES UNI femoral component cannot be combined with CERAVER cemented patellar buttons or extensions.
- B. The tibial component consists of two parts (a baseplate and an insert) which must be assembled at the time of the surgical procedure; **The engagement mechanism can only be active once. Do not remove and reinsert the polyethylene inserts once they are in place.**
- C. CERAVER components, instruments and test prostheses should not be used with those of another manufacturer.

## 7 - PRECAUTIONS

### 7.1 Preoperative precautions

- A. When implanting a unicompartamental knee prosthesis is being considered, particularly in a young and/or active person, all aspects of the surgery and the implant should be discussed with the patient before the procedure. This consultation should cover the limitations of joint reconstruction, any patient-specific limitations, the possible consequences of these limitations and thus the absolute requirement to follow the doctor's preoperative instructions.
- B. During the preoperative period, it is necessary to consider and eliminate any risks of allergy and other reactions to the implant materials, however rare these risks may be.
- C. Surgical technique and preoperative planning for the implantation of a HERMES UNI unicompartamental knee prosthesis correspond to the basic principles of good surgical practices for the fitting of a unicompartamental knee prosthesis: an in-depth knowledge of the surgical technique is required.
- D. Specific surgical tools must be used and the surgeon must be familiar with them. Cutting alignment and guiding instruments should be checked before surgery. Damaged or bent instruments may lead to incorrect implant positioning, thus compromising the success of the procedure. The technical surgical brochure associated with the type of prosthesis is supplied with each set of instruments and is available from CERAVER.

• **HERMES UNI** - Instructions for use

---

- E. Radiographic templates should be used to estimate the size and positioning of implants as well as its alignment to the HKA axis. A sufficient stock of implants of various sizes must be available at the time of surgery, including sizes smaller and larger than planned. The integrity of all packaging and implants must be carefully checked before surgery.
- F. It is up to the surgeon to gage the safety and efficacy of bilateral HERMES UNI prosthesis implantation.
- G. **Prior to each operation, it is necessary to ensure that a more constrained prosthesis is available.**

## 7.2 Perioperative precautions

- A. If a contraindication comes to light during the operation, it is up to the surgeon to decide whether or not to continue with implantation or implant another type of prosthesis.
- B. The tibial component should be implanted so as to ensure adequate coverage of the bone's cortical substance in all directions (medial-lateral and anterior-posterior).
- C. It is important to ensure adequate tibial surface preparation and cleaning to achieve correct prosthesis anchoring. Keep bone ablation to the strict minimum to allow implant placement. Excessive bone ablation can lead to mechanical disorders, bone resorption and, consequently, failure of the procedure as a result of loosening or deformation of the implant. When preparing and placing the tibial components, correct tibial alignment should be checked.
- D. The use of polymethyl methacrylate (PMMA) cement helps fix, support and stabilize certain bone implants but it cannot either replace the support function of healthy bone tissue or preclude the need for additional support during the consolidation period. When cement is used, complete support of the cement on all parts of the component included in it should be checked in order to prevent the concentration of stresses that could compromise the success of the procedure.
- E. Ablation of osteophytes is also an important parameter that may affect the implant's kinematics and particularly the degree of flexion.
- F. Before closing the surgical site, the knee should be cleaned carefully to remove all bone fragments, ectopic bone tissue, cement, etc. The presence of foreign particles in the metal/plastic joint interface may lead to excessive wear and/or friction. Ectopic bone tissue and/or bone spurs may cause dislocation or painful and restricted mobilization. Check joint mobility carefully to detect any early contact or instability.

## 7.3 Postoperative precautions

- A. The patient must be told about the restrictions associated with a HERMES UNI unicompartmental knee prosthesis and warned about using the operated limb for support and informed of permitted activity levels. Monitor initial support attempts closely.
- B. Carefully structure initial postoperative care in such a manner as to maintain joint mobility and to prevent any dislocation or thromboembolism.
- C. Structure postoperative treatment, patient care (dressing changes, bed pan use, etc.) and physical activities such as to prevent any overload of the operated knee. The patient's age and/or bone quality may require the limited load period to be extended.
- D. All patients implanted with a CERAVER unicompartmental knee prosthesis should see an orthopedic surgeon for postoperative follow-up. This examination is important to detect any signs of abnormal prosthetic wear or loosening before the appearance of functional manifestations.
- E. Encourage the patient to notify his/her doctor immediately of any unusual changes to the operated joint.

- F. When annual examination is impossible, the surgeon should receive, at least once a year, a control radiograph for evaluation.
- G. The implants are compatible with MRI examination. They may nonetheless create artefacts that should be taken into account for a perfect interpretation of the examination. However, a period of 48 hours is recommended between implant placement and MRI examination to avoid any movement of the device. Safety, compatibility, heating and migration of HERMES UNI implants have not been tested in an MRI. An assessment has been performed in comparison to similar devices available on the market.

## 8 - ADVERSE EFFECTS

The most commonly reported potential complications are:

- A. Anesthesia complications;
- B. Bone cement implantation syndrome;
- C. Regressive or permanent vascular, musculoligamentous or nerve lesions;
- D. Per and postoperative fractures;
- E. Hemorrhage, heterologous/autologous transfusion;
- F. Misalignment of prosthetic components;
- G. Other complications associated with a general surgical procedure (digestive, cardiovascular, pulmonary, urological problems, etc.), medicines or ancillary surgical equipment used, blood, etc...
- H. Hematomas, sanguineous or non-sanguineous intra-articular effusions;
- I. Thromboembolic complications: phlebitis and pulmonary embolism;
- J. Wound healing disorders;
- K. Early and late-onset infections;
- L. Inflammation;
- M. Algodystrophy;
- N. Ectopic ossification;
- O. Persistent knee pain;
- P. Reduced or low amplitude of movement;
- Q. Knee joint instability, femoropatellar or femorotibial dislocation or subluxation;
- R. One or more torn cruciate ligaments requiring a new operation;
- S. Allergy to implant components;
- T. Spread of arthrosis to the femorotibial compartment without the prosthesis and/or the femoropatellar compartment;
- U. Reaction to products of corrosion, free particles of cement or wear that may cause osteolysis;
- V. Movement, loosening or breakage of prosthesis components;
- W. New operation.

## 9 - PACKAGING AND STORAGE

- A. The prostheses are stored in their original packaging.
- B. The various components are single-use, they are unit packaged in double bags and provided sterile.
- C. Do not use the device after the expiry date. Contact your CERAVER representative to return the product.

• **HERMES UNI** - Instructions for use

---

**10 - STERILIZATION**

- A. All CERAVER knee implants are sterilized using gamma radiation at a minimum dose of 25 kGy; UHMWPE components are argon-conditioned to avoid polymer oxidation during sterilization.
- B. Before use, it is important to check the integrity of the sterility protection.

**11 - RESTERILIZATION AND REUSE**

- A. Resterilization and/or reuse of any of the components of a CERAVER knee replacement is strictly forbidden because of the risk of contamination that could lead to infection.
- B. Resterilization and/or reuse of any of the components of a CERAVER knee replacement may compromise the device's mechanical performance and biological safety. Never reuse surgical implants. If an implant is removed after it has been implanted, it must be disposed of. Although it may appear intact, it may exhibit minor faults and signs of internal stresses that may undermine its integrity. Only new implants should be used and their configuration must not be modified before use.

**12 - DISPOSAL AFTER USE**

If a device is removed for follow-up replacement, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

**13 - NOTE**

The HERMES UNI unicompartmental femorotibial prosthesis was granted CE Marking in 1990.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

**CERAVER® - Les Laboratoires Osteal Médical**  
69 rue de la Belle Etoile  
CS. 54263 - 95957 Roissy CDG Cedex - France  
Tél. : 33 (0)1 48 63 88 63 - Fax : 33 (0)1 48 63 88 99  
[WWW.CERAVER.COM](http://WWW.CERAVER.COM)

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