# **BOUTON ROTULIEN CIMENTÉ**

STERILE R

CE 0459

- **BOUTON ROTULIEN CIMENTÉ CERAVÉR**  
Notice d'instructions
- **CEMENTED CERAVÉR PATELLAR COMPONENT**  
Instructions for use
- **ZEMENTIERTE PATELLAKOMPONENTE CERAVÉR**  
Anleitung
- **COMPONENTE ROTULIANO CEMENTADO CERAVÉR**  
Manual de instrucciones
- **BOTTONE ROTULEO CEMENTATO CERAVÉR**  
Istruzioni d'uso
- **НАДКОЛЕННИК ЦЕМЕНТНОЙ  
ФИКСАЦИИ CERAVÉR**  
Инструкция по эксплуатации

LES LABORATOIRES OSTEAL MEDICAL

STERILE R

## • BOUTON ROTULIEN CIMENTÉ CERAVÉR

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STERILE R

0459

**• CEMENTED CERAVER PATELLAR COMPONENT**

Instructions for use

Instructions for using the following component:

Cemented CERAVER patellar component	5090	5091	5092	5093
-------------------------------------	------	------	------	------

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken
	Batch code		
	Sterilized by irradiation		CE marking followed by the code for the notified body (0459, LNE/G-MED)
	Catalogue reference		

**Before using a CERAVER patellar component, the attention of the surgeon is drawn to the following recommendations:**

### **1 - GENERAL POINTS**

- A. The patellar component is an optional component of a CERAVER knee prosthesis and cannot be used on its own. Ceraver recommends fitting the patella with a prosthesis. However, according to the characteristics of the patient (bone quality, age, physical activity, weight), it is up to the surgeon to decide not to implant the prosthesis.  
The use of the CERAVER patellar component with components of other origins is strictly prohibited.  
The manufacturer and supplier of a CERAVER patellar component guarantee the quality of the manufacturing process and material used. A CERAVER patellar component already used must not be implanted again.
- B. Several criteria are key to the success of a procedure to implant a CERAVER knee prosthesis:
- Appropriate patient selection for the procedure;
  - Correct implant selection;
  - Observance of the operative technique and correct use of the surgical instruments specially designed by CERAVER;
  - Application of strict aseptic technique.
- C. The surgeon is responsible for any complications that may arise from erroneous indication and/or use of the instruments associated with the operative technique and/or of poor aseptic practice; neither the manufacturer nor the supplier of the CERAVER products can be blamed for these complications.
- D. The surgeon must be aware of all problems posed by implanting a knee prosthesis. The surgical procedure for implantation of a CERAVER total knee prosthesis can be best learned at a clinic or training centre with surgeons who have acquired expertise in the technique.
- E. The surgeon must inform patients receiving prosthetic knee joints of the fact that their weight and activity levels may influence the durability and performance of the prosthesis.**

### **2 - DESCRIPTION OF THE CERAVER PATELLAR COMPONENT**

The CERAVER patellar component (UHMWPE as per NF ISO 5834-2) is a cemented knee prosthesis component. The centre to centre distance between the anchoring points is identical, whatever the size of the patellar component.

Reference	Diameter (mm)	Thickness (mm)
5090	Ø31	8
5091	Ø34	9
5092	Ø36	10
5093	Ø39	11

### 3 - INDICATIONS

Refer to the instructions for use of the implant associated with the patellar component.

### 4 - ANTICIPATED PERFORMANCE

Refer to the instructions for use of the implant associated with the patellar component.

### 5 - CONTRAINDICATIONS

Refer to the instructions for use of the implant associated with the patellar component.

### 6 - WARNINGS

Refer to the instructions for use of the implant associated with the patellar component.

The following factors can be extremely important for the success of the operation:

#### A. Implant selection

Careful selection of the right implant size is necessary for correct joint function. The degree of success of a total knee prosthesis implantation depends on choosing an implant of the right size and model. The patient's age, level of physical activity, weight, muscle condition and bone quality must also be considered.

#### B. Handling the medical device

The exposed implant joint surfaces must not be placed in contact with any metallic objects or other hard objects that could potentially damage the surface.

#### C. Assembly and compatibility of components

The patellar component is fully compatible with all the femoral components of the ranges:

- Cemented HERMES PS Revision;
- Cemented HiFit PS;
- Cemented HERMES FP unicompartmental femoropatellar prosthesis.

### 7 - PRECAUTIONS

Refer to the instructions for use of the implant associated with the patellar component.

### 8 - ADVERSE EFFECTS

Refer to the instructions for use of the implant associated with the patellar component.

### 9 - PACKAGING AND STORAGE

- A. The patellar components are stored in their original packaging.
- B. The various components are single-use, they are unit packaged in double bags and provided sterile.
- C. Do not use the device after the expiry date. Contact your CERAVÉR representative to return the product.

## 10 - STERILIZATION

- A. All CERAVER knee implants are sterilized using gamma radiation at a minimum dose of 25 kGy; UHMWPE components are argon-conditioned to avoid polymer oxidation during sterilization.
- B. Before use, it is important to check the integrity of the sterility protection.

## 11 - RESTERILIZATION AND REUSE

- A. Resterilization and/or reuse of any of the components of a CERAVER knee replacement is strictly forbidden because of the risk of contamination that could lead to infection.
- B. Resterilization and/or reuse of any of the components of a CERAVER knee replacement may compromise the device's mechanical performance and biological safety. **Never reuse surgical implants.** If an implant is removed after it has been implanted, it must be disposed of. Although it may appear intact, it may exhibit minor faults and signs of internal stresses that may undermine its integrity. Only new implants should be used and their configuration must not be modified before use.

## 12 - DISPOSAL AFTER USE

If a device is removed for revision surgery, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

## 13 - NOTE

The Ceraver knee prosthesis patellar component was granted CE Marking in 1996.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

**CERAVER - Les Laboratoires Osteal Médical**  
69 rue de la Belle Etoile  
CS. 54263 - 95957 Roissy CDG Cedex - France  
Tel.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99  
[WWW.CERAVER.COM](http://WWW.CERAVER.COM) ; [adv@ceraver.com](mailto:adv@ceraver.com)

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