

The bone substitute for dental and orthopaedic applications

CALCIRESORB

Tricalcium phosphate β

Classic

- 1991 : the first implantation
- To date : hundreds of thousands implantations



Dental Bone substitutes
Fuse, Rebuild, Strengthen

 **CERAVER**
EXPERIENCE - INNOVATION

CALCIRESORB

Classic

Bone substitute material is intended to be used in maxillofacial and dental surgery for :

- Bone defects or cavities requiring filling
- Alveolar filling post-extractional
- Elevation of the sinus floor

Calciresorb helps you to optimize surgical procedures to foster success of your future implants with :

- **Full and progressive resorption** : renews bone integrity in 3 to 6 months
- **Natural regeneration of neo-formed bone** : micro and macroporous structure with osteo conductive properties
- **Restoration of bone volume** : calcium and phosphate ions liberation helping new bone growth
- **Adapted volume** : even for extensive lesions – 2cc fragments
- **Ease of implementation** : through simple blood or physiological liquid suction – 0.8 cc syringe

COMPOSITION	SIZE	POROSITY	PORE SIZES	RESORPTION	REF
β tricalcium phosphate	Fragments of 500 μm to 1 mm	65 %	100 < Ø _{macropores} < 400μm	3 to 6 months according to patient's physiology	4633 0,8 cc fragments pre-filled in syringe
					4635 2 cc fragments in vial

0,8 cc fragments pre-filled in syringe



2 cc fragments in vial



Document intended for the exclusive use of healthcare professionals. CALCIRESORB® CLASSIC – bone substitutes – is a class III CE marked medical device made by CERAVÉR – LES LABORATOIRES OSTEAL MEDICAL Company and for which Conformity assessment was carried out by Notified Body G-MED n°0459. CALCIRESORB is intended to bone-void filling by vascular invasion of the ceramic pores and by progressive osseointegration which leads to the replacement of the implant by the bone. Before any surgical procedure, read carefully instructions for use. CALCIRESORB® CLASSIC is a medical device covered by Health Insurance under certain conditions. Consult the terms and conditions on ameli.fr website.

CERAVÉR
EXPERIENCE - INNOVATION

CERAVÉR – LES LABORATOIRES OSTEAL MEDICAL
69, RUE DE LA BELLE ÉTOILE
95957 ROISSY CDG CEDEX – FRANCE
TEL. : (+33) 1 48 63 88 63 - FAX : (+33) 1 49 89 05 05
WWW.CERAVÉR.COM – CONTACT@CERAVÉR.COM

CERAVÉR OSTEAL AMERICA INC.
2277 GUENETTE
SAINT LAURENT QUEBEC
H4R 2E9 CANADA
TEL: 19 1 514 856 24 35
FAX: 19 1 514 856 24 44
coai@qc.aira.com

CERAVÉR IBERICA
BELTRAN BAGUENA, 5-10-10
46009 VALENCIA
ESPAÑA
TEL: +34 96 348 10 17
FAX: +34 96 348 10 02
info@ceraver.es

CERAVÉR ITALIA
VIA PRATI, 29
40012 CALDERARA DI RENO (BO).
ITALIA
TEL: +39 051 72 88 76
FAX: +39 051 72 78 03
info@ceraver.it

CERAVÉR SOUTH AFRICA (PTY) LTD.
11 POLO CRESCENT,
WOODMEAD OFFICE PARK,
WOODMEAD, 2191, SANDTON
GAUTENG - SOUTH AFRICA
TEL: +27-11-656-3719
FAX: +27-865-474-635

