


OSTEAL



Les Tiges de Reprise OSTEAL The OSTEAL Revision Stem


 **CERAVER**
EXPERIENCE - INNOVATION

GAMME - RANGE

Les tiges longues OSTEAL complètent la gamme standard pour répondre à certains cas de reprise.

OSTEAL long stems complete the standard range for cases of revision.

LARO STANDARD
(150 mm)



LARO DE REPRISE
(250 mm)



PROFIL PROFILE	LONGUEUR (MM) LENGHT (MM)	RÉFÉRENCES REFERENCES	MATÉRIAUX MATERIAL	NORMES STANDARDS
DND	200 mm	10 004	Alliage de Titane TiAl6V4 Titanium alloy TiAl6V4	ISO5832-3
	250 mm	10 005		
LAR	200 mm	10 007		ASTM F136
	250 mm	10 008		
	300 mm	10 009		

Document destiné exclusivement aux professionnels de Santé. Tige OSTEAL® - tige fémorale cimentée - est un dispositif médical marqué CE de classe III fabriqué par la société CERAVER – LES LABORATOIRES OSTEAL MEDICAL. L'évaluation de la conformité a été réalisée par l'Organisme Notifié G-MED N°0459. La tige OSTEAL est destinée à remplacer partiellement ou totalement une articulation de hanche ne pouvant être soignées par le biais d'autres thérapies. Avant toute opération chirurgicale, lisez attentivement la notice d'instructions, et la technique opératoire. Pour le bon usage et la pose de ces dispositifs, les professionnels qualifiés doivent utiliser les instruments du kit associé. Tige OSTEAL® est un dispositif médical pris en charge par l'Assurance Maladie dans certaines conditions. Consulter les modalités sur le site ameli.fr.

Document intended for the exclusive use of healthcare professionals. OSTEAL Stem® - cemented femoral stem – is a class III CE marked medical device made by CERAVER – LES LABORATOIRES OSTEAL MEDICAL Company and for which Conformity assessment was carried out by Notified Body G-MED n°0459. The OSTEAL stem is intended to replace completely a hip joint that cannot be treated through other therapies. Before any surgical procedure, read carefully instructions for use and surgical technique. For proper use and installation of these devices, qualified professionals must use instruments of the associated kit. is a medical device covered by Health Insurance under certain conditions. Consult the terms and conditions on ameli.fr website.



CERAVER – LES LABORATOIRES OSTEAL MEDICAL
69, RUE DE LA BELLE ETOILE
95957 ROISSY CDG CEDEX – FRANCE
TEL. : (+33) 1 48 63 88 63 - FAX : (+33) 1 49 89 05 05
WWW.CERAVER.COM – CONTACT@CERAVER.COM

CERAVER OSTEAL AMERICA INC.
2277 GUENETTE
SAINT LAURENT QUEBEC
H4R 2E9 CANADA
TEL: 19 1 514 856 24 35
FAX: 19 1 514 856 24 44
coai@qc.aira.com

CERAVER IBERICA
BELTRAN BAGUENA, 5-10-10
46009 VALENCIA
ESPAÑA
TEL: +34 96 348 10 17
FAX: +34 96 348 10 02
info@ceraver.es

CERAVER ITALIA
VIA PRATI, 29
40012 CALDERARA DI RENO (BO).
ITALIA
TEL: +39 051 72 88 76
FAX: +39 051 72 78 03
info@ceraver.it

CERAVER SOUTH AFRICA (PTY) LTD.
11 POLO CRESCENT,
WOODMEAD OFFICE PARK,
WOODMEAD, 2191, SANDTON
GAUTENG - SOUTH AFRICA
TEL: +27-11-656-3719
FAX: +27-865-474-635

